

Egyptian Society of
CARDIOLOGY

Mansoura University Tanta University Zagazig University



45TH Annual International Congress of the
EGYPTIAN SOCIETY OF CARDIOLOGY
CardioEgypt 2018


Coarctation angioplasty and stenting step by step



Maiy H. El Sayed, MD
Professor and Head of Cardiology Department
Ain Shams University
Cairo, Egypt



- Coarctation of Aorta accounts for approximately 7% of all congenital cardiac defects, with high occurrence in males
- Associated arch hypoplasia is common
- Associated lesions : bicuspid aortic valve, PDA, VSD, anomalous right subclavian artery




45TH Annual International Congress of the
EGYPTIAN SOCIETY OF CARDIOLOGY
CardioEgypt 2018

www.cardioegypt.com



Diagnosis and patient selection

- Clinical examination
- Echocardiography
- MSCT
- MRI

Echocardiography



MSCT



45th Annual International Congress of the
EGYPTIAN SOCIETY OF CARDIOLOGY
CardioEgypt 2018

Indications of Balloon Angioplasty/ stenting

Balloon Angioplasty

- Neonatal period and/ Or first 6 months of life in case of poor LV systolic function, sepsis, hemorrhage
- Discrete coarctation whether native or S/P surgery < 15 Kg
- Children below 6-8 years
- Adults with favorable anatomy.

Stenting

- Native or s/p surgery - weight > 15-20 Kg
- Long segment coarctation
- Associated hypoplasia of isthmus and/or aortic arch
- Tortuous coarctation
- Recurrent after surgery or aneurysm following prior balloon or surgery



45th Annual International Congress of the
EGYPTIAN SOCIETY OF CARDIOLOGY
CardioEgypt 2018

www.cardioegypt.com

Equipment – neonates and children

- 4-5F introducer
- Multipurpose, pigtail, JR
- 0.018, 0.025, 0.035 GW
- Tyshak mini in neonatal period
- Older children – Tyshak II, Opta, Z-med, Powerflex



Tyshak mini Double Tapered Balloon

- Super thin for a low deflated profile that maintains tip flexibility.
- Exceptionally low profile
- Requires the smallest introducer possible.



Tyshak, TyshakII


- Coaxially constructed catheter with a distally mounted non-compliant balloon.
- Extremely low profile.
- Due to the coaxial construction, the catheter inflates and deflates extremely fast.

Stents for neonates and small children



Formula Balloon expandable stent Cook


- Open cell
- Low crossing profile
- No shortening
- 5-10mm can be redilated
- 4-5F flexor sheath – 6F coronary guide cath



VALEO Balloon expandable stent BARD

- Uniquely flexible stent
- Competitive radial strength with minimal recoil
- 6-10mm can be redilated - forshortens
- 6-7F sheath – 7F coronary guide cath

Equipment - Stenting



The **Z-MED™** high pressure dilatation catheter is a coaxially constructed catheter with a distally mounted non-compliant high pressure balloon.

0.035" stiff wire




BIB balloon

An inner balloon 1/2 of the balloon diameter of the outer balloon, 1cm shorter and has a burst pressure of 45.5atm.

Equipment

Palmaz and Genesis Stents



- Balloon expandable stainless steel.
- Closed cell design which gives them high radial strength but makes them less flexible.
- Hand mounted and crimped on the delivery balloon.

Equipment

CP Stent™




- Composed of 0.013" platinum/iridium wire that is arranged in a "zig" pattern, laser welded at each joint and then over brazed with 24K gold.
- It allows expansion from 12.0mm to 24.0mm.

The **Covered CP Stent™** is comprised of the Bare CP Stent that is covered with an expandable sleeve of ePTFE.

Hand crimped or balloon mounted

Balloon Angioplasty – step by step


- Femoral artery access – 18G needle or 19G scalp
 - Anatomic or ultrasound guided
- <10kgs 3-4F sheath
- Heparin 100 units/kg
- Crossing Multipurpose catheter or JR
- PTCA wires (hydrophilic)
- Pressure gradient and Angiography – pigtail catheter
- Balloon positioning – inflation
- Final Pressure gradient and angiogram
- Hemostasis to avoid complications



Femoral artery access

- The bifurcation of the CFA occurs in zones 2, 3, 4 and 5, 1%, 9%, 43% and 47% of the time, respectively
- Occurs within the lower third of the femoral head or below the lower border of the femoral head in 90% of patients

A 18 G needle could be ideal
Puncture the skin at zone 5 inferior border of head of femur. Enter the artery at mid point in the Zone 3.
The chances of hitting the femoral artery is near 95 %




Cardiovascular Intervention and Therapeutics January 2014, Volume 29, 18-23 Madjid Chinikar

Step by Step Stenting

- Femoral artery access – heparin 100U/KG – ACT > 200
- Angiogram in descending aorta
- Crossing the coarctation
- Angiogram above coarctation – measurements
- Crossing with the stiff wire
- Long mullins sheath [size depending on balloon and stent used]
- Preparing balloon and crimping the stent
- Confirming position
- Stent deployment
- Withdrawing the balloon inside sheath
- Final pressure gradient and angiogram

Size of balloon and stent



CardioEgypt 2018

www.cardioegypt.com

Long sheath



45TH
45th Annual Meeting
EGYPTIAN SOCIETY OF CARDIOLOGY
CardioEgypt 2018

www.cardioegypt.com

Preparing balloon




Crimping the stent



45th
EG
CardioEgypt 2018

www.cardioegypt.com

Stent positioning - Rotational angio



EGYPTIAN SOCIETY OF CARDIOLOGY
CardioEgypt 2018

In conclusion

- Balloon angioplasty in children, and stent deployment in adolescents and adults are considered the standard of care
- The anatomy of the coarcted aortic segment and the aortic arch would greatly influence the method of therapy in a given patient.
- Depending on the experience of the center most coarctations can be treated in the cath lab with only cases with very complex anatomy to be referred to surgery

45TH
45TH Annual International Congress of the
EGYPTIAN SOCIETY OF CARDIOLOGY
CardioEgypt 2018

www.cardioegypt.com

CSI Africa 2018 Cairo, Egypt

November 30 – December 1, 2018

For more information please visit the website:

<http://www.csi-congress.org/csi-africa.php>

